

BANDAGING AND SUTURE



STUMP BANDAGING SIMULATORS

Practice and demonstrate standard bandaging procedures. The lower torso is long enough to allow carrying the bandage around the body over the uninvolved hip at the level of the iliac crest. (NOTE: There is no embedded, palpable skeletal material.) Both the upper and lower portions can be used to demonstrate the attachment of prosthetic devices. The upper torso includes the two arms — one is amputated above the elbow; the other above the wrist. Both arms are slightly extended in a patient-like position to facilitate bandaging. The compressibility of the material very closely duplicates actual experience with a patient. The lower torso is representative of a patient in a supine position with legs slightly abducted. One leg is amputated below the knee; the other at mid-thigh.

SBM-39

Upper Stump Bandaging Simulator

Includes an instruction booklet. Individually boxed.
Dimensions: 66 x 46 x 33 cm.
Sh. wt. 7.50 kgs.

SBM-43


Lower Stump Bandaging Simulator

Includes an instruction booklet. Individually boxed.
Dimensions: 89 x 38 x 25 cm.
Sh. wt. 6.80 kgs.

SBM-46

Set of 2 Stump Bandaging Simulators

Includes instruction booklets. Individually boxed.
Dimensions: 89 x 38 x 25 cm.
Sh. wt. 11.80 kgs.



SVQ-30

"Surgical Sally" Bandaging Simulator

A readily available patient with 14 surgical wounds including a Mid-Sternal Split with 2 simulated drains, a Sacral Decubitus Ulcer - Stage 2, and a Leg Amputation Stump. The skin has also been reformulated so it is drier, allowing the bandages to adhere better. Wound closures feature surgical staples for the Mid-Sternal Split, Thoracotomy, Nephrectomy, Laparotomy, Abdominal Hysterectomy, and Amputation Stump. Staples and sutures cannot be removed. Surgical Sally is a replica of a 29" x 17" female. Special vinyl formulations have been used to create the most lifelike look and feel possible. The flexible, flesh-colored skin realistically responds to adhesives and all types of bandaging procedures. Three year warranty. Surgical Sally features the following simulated wounds:

- Thyroidectomy • Mid-Sternal Split - With Chest Tube Drains
 - Mastectomy - With Simulated Drain
 - Cholecystectomy - With Simulated T-Tube
 - Laparotomy • Appendectomy • Colostomy • Ileostomy
 - Abdominal Hysterectomy
 - Thoracotomy • Nephrectomy • Laminectomy
 - Sacral Decubitus Ulcer - Stage 2 • Leg Amputation Stump
- The Surgical Bandaging Simulator comes complete with hard carrying case.
Dimensions: 83 x 63 x 30 cm.
Sh. wt. 11.80 kgs.



Suture And Stapling Practice Trainers

Made with a soft vinyl skin over a core of stitchable foam to provide a lifelike suturing experience for students or a realistic suturing demonstration by an instructor. Realistic skin texture with wrinkles, pores, and visible fingerprints. Soft and pliable for easy sewing, the skin is tough enough that sutures will not pull out when tightened. The unit is provided with three "wounds." These "wounds" can be sutured repeatedly until the skin around them is finally worn out; then, new "wounds" can be cut. Although disposable, the arm provides hundreds of suturing experiences before wearing out. Conservatively, over one hundred cuts can be made on the trainer and each of these cuts can be sutured several times. The Suture Practice Leg is complete with a starter suturing kit and instruction booklet.

BPS-45

Suture And Stapling Practice Arm

Complete with a starter package of sutures and instruction booklet.
Sh. wt. 1.10 kgs.
Dimensions: 64 x 18 x 13 cm.

PPS-55

Suture And Stapling Practice Leg

Complete with a starter package of sutures and instruction booklet.
Sh. wt. 2.00 kgs.

